

DIGITAL LAB NOTEBOOK OF KEVIN MURRAY

HONOUS PROJECT, 2013

Jointly supervised by Justin Borevitz and Barry Pogson

Last updated 18:27

2013-01-21

1 Practice RNA Extraction

1.1 Quantification of RNA samples

Aim

- Determine qty of RNA in previously extracted samples

Method

- Nanodropped RNA extraction from 15/1/13??
- Standard protocol, used sterile milliQ water as blank.

Result

- Of the 14 samples, 10 had reasonable amounts of RNA, and 260/280 ratios were above 1.8 in all but one case. (see ./jan/20130121-PracticeRNASamples.ods)

Attachments

- ./jan/20130121-PracticeRNAExtractionSamples.csv
- ./jan/20130121-PracticeRNAExtractionSamples.ndv
- ./jan/20130121-PracticeRNASamples.ods