

DIGITAL LAB NOTEBOOK OF KEVIN MURRAY

HONOUS PROJECT, 2013

Jointly supervised by Justin Borevitz and Barry Pogson

Last updated at 16:55 on Wednesday 30th January, 2013

Mon 2012-12-10

1 Final Harvest of Keng's RIX lines

1.1 Aim

Harvest lines after 1 week of HL stress.

1.2 Method

- An Eppendorf 1.2mL deep well plate was placed on dry ice for ≈ 10 minutes before sampling to allow to cool.
- Whole leaves were excised and placed into 1.2mL Eppendorf 96 deep well plate.
- Where possible, the largest mature leaf was taken. In some cases, this was hard to determine, so the youngest of the fully-expanded leaves was taken (as this was generally also the largest leaf). Some plants were very small, and had only juvenile leaves, in which case the largest juvenile leaf was taken.

1.3 Results

The following table describes the plate layout.

		Well	Line	Comments	Well	Line	Comments
A1	87	2nd plant with orange label			A7	65	NPQ4
B1	100				B7	71	
C1	92				C7	53	
D1	OSB2				D7	55	
E1	93				E7	56	
F1	94				F7	50	
G1	1*cvi				G7	51	
H1	99				H7	52	
A2	63				A8	78	
B2	98				B8	80	
C2	98				C8	85	
D2	70				D8	74	
E2	67				E8	73	
F2	99				F8	76	
G2	66				G8	83	
H2	100				H8	81	
A3	65				A9	89	
B3	72				B9	90	
C3	12				C9	91	
D3	11				D9	16	
E3	10				E9	17	
F3	12				F9	18	

Well			Line			Comments		
G3	9					G9	21	
H3	8					H9	23	
A4	7					A10	25	
B4	6					B10	57	
C4	5					C10	28	
D4	13					D10	29	
E4	2					E10	30	
F4	1*cvi					F10	31	
G4	49					G10		
H4	47					H10		
A5	46					A11		
B5	42					B11		
C5	45					C11		
D5	71					D11		
E5	39					E11		
F5	40					F11		
G5	41					G11		
H5	43					H11		
A6	38					A12		
B6	39					B12		
C6	36					C12		
D6	33					D12		
E6	61					E12		
F6	62					F12		
G6	68					G12		
H6	58					H12		

Attachments:

- dec12/20121210-harvest_photos.tar.bz2
- MD5SUM:40dae2cad3babaa3c32f0d35a9d9442c

Mon 2013-01-21

1 Quantification of RNA samples

1.1 Aim

- Determine qty of RNA in previously extracted samples

1.2 Method

- Nanodropped RNA extraction from 15/1/13??
- Standard protocol, used sterile milliQ water as blank.

1.3 Result

- Of the 14 samples, 10 had reasonable amounts of RNA, and 260/280 ratios were above 1.8 in all but one case. (see ./jan/20130121-PracticeRNASamples.ods)

1.4 Attachements

- ./jan/20130121-PracticeRNAExtractionSamples.csv
- ./jan/20130121-PracticeRNAExtractionSamples.ndv
- ./jan/20130121-PracticeRNASamples.ods

2 MADE: 10x MOPS Solution

Method

- Add 41.8g RNA only MOPS to beaker
- Add 450mL DEPC H₂O, mix w/ stirrer bar on mag stirrer
- Add 26.6mL 3M Sodium Acetate (0.22um Filtered before use)
- Add 10mL RNA only 0.5M EDTA
- pH to 7 with 5M NaOH
- Top up to 500 mL with DEPC H₂O
- Use 10ml per 100mL MOPS gel

Mon 2013-01-21

3 MADE: RNA Denaturing Gel (MOPS)

Method

- Melt 1g RNase-free Agarose in 72ml DEPC H₂O
- Add 10mL 10x MOPS
- Add 18mL 37% Formaldehyde
- Pour in RNA-only gel tank, previously washed with 0.5% SDS and RNase-zap

Tue 2013-01-22

1 Denature RNA for RNA gels

1.1 Method

- Dilute RNA to 100ng/uL
- Add RNA gel loading buffer (Obtained from Pete Crisp)
- Incubate at 65 degrees for 10 minutes. The samples were incubated for 10 minutes on the evening of 2013-01-21, but the gels were not run until 2013-01-22, so they were denatured for a further 2 minutes at 65 degrees

2 TBE Gel

2.1 Aim

- To compare TBE and denaturing/MOPS gels for RNA

2.2 Method

- Dissolve 1g RNAase-free agarose in 90mL DEPC water
- Add 10mL RNAase-free TBE (prepared using DEPC Water, obtained from Pete Crisp)
- Pour in RNA-only gel tank, previously washed with 0.5% SDS or RNAse-zap
- Then, load denatured samples, and run in RNAse-free 1x TBE
- Run at $\approx 80V$, $\approx 40-50mA$ for $\approx 1.75h$
- Stain gel in 0.5ug/ml Ethidium Br in DEPC water?? for 10 min on orbital shaker, and photograph.

2.3 Result

See Figure ?? below.


Gel indicates some degradation of RNA, however most samples are OK. Sample order is (left to right) A2, A3, A5, A6, A7, B3, B5, B7. A7 appears to have no RNA, although this is probably a misloading error. Overall, the TBE gel appears to be of more use than the MOPS gel.

3 MOPS gel

3.1 Aim

- Determine quality of RNA and Compare MOPS with TBE for RNA gels

Figure 1: TBE Gel of Practice RNA samples, 2013-01-22




3.2 Method

- Load samples after denaturing as above. Sample order is (left to right) A2, A3, A5, A6, A7, B3, B5, B7.
- Run gel in RNase free 1x MOPS at $\approx 80V$, $\approx 100mA$ for $\approx 1.75h$ as per TBE gel above.
- Stain gel in $0.5\mu g/ml$ Ethidium Br in DEPC water?? for 10 min on orbital shaker.
- Destain on orbital shaker in 1x MOPS, and photograph. Gel disintegrated whilst destaining.

3.3 Results


Figure 2: MOPS Gel of Practice RNA samples, 2013-01-22



See Figures ?? and ??

Mops gel confirms that the rna was of reasonable quality. The MOPS gel appears to be of less use than the TBE gel.

Figure 3: MOPS Gel of Practice RNA samples, 2013-01-22. Contrast adjusted.



Tue 2013-01-29

1 Seed Stock Levels

The stocks of Joost's RIX set were checked. Seed lines were classified as having either plenty (+), limited(?) or no (-) seed. The levels of each line are shown in the table below.

	Line	Desc	Count	Line	Desc	Count
1	Col-0 ₁ 915	not in box		51	168 x 22	+
2	Col-0 ₄ 936	not in box		52	169 x 175	not in box
3	Cvi x Cvi	-		53	17 x 21	+
4	1 x Cvi	not in box		54	170 x 24	-
5	1 x 146	+		55	171 x 143	+
6	10 x 26	+		56	174 x 34	+
7	101 x 176	+		57	180 x 157	?
8	Ll-1	not in box		58	183 x 118	-
9	105 x 145	not in box		59	186 x 27	+
10	107 x 124	not in box		60	187 x 190	-
11	109 x 185	not in box		61	187 x 69	not in box
12	109 x 47	not in box		62	189 x 133	-
13	110 x 32	?		63	19 x 173	+
14	112 x 30	+		64	19 x 67	+
15	113 x 141	-		65	190 x 176	+
16	114 x 3	?		66	191 x 31	not in box
17	114 x 60	+		67	192 x 189	+
18	115 x 126	+		68	20 x 138	+
19	117 x 73	?		69	21 x 22	-
20	118 x 108	+		70	24 x 171	+
21	118 x 164	+		71	25 x 9	+
22	119 x 177	-		72	26 x 74	+
23	12 x 142	+		73	33 x 58	not in box
24	122 x 42	?		74	35 x 120	-
25	125 x 117	+		75	38 x 35	+
26	128 x 6	+		76	39 x 27	not in box
27	132 x 129	not in box		77	40 x 74	-
28	133 x 35	+		78	npq4	not in box
29	134 x 29	+		79	43 x 131	not in box
30	135 x 10	+		80	44 x 50	+
31	135 x 140	?		81	45 x 23	+
32	136 x 102	+		82	46 x 29	+
33	165 x 137	not in box		83	48 x 160	+
34	139 x 162	-		84	49 x 158	+
35	139 x 36	+		85	5 x 172	not in box
36	14 x 4	+		86	5 x 188	not in box
37	146 x 64	not in box		87	51 x 111	+
38	147 x 50	+		88	51 x 18	+
39	147 x 69	+		89	54 x 183	+

	Line	Desc	Count	Line	Desc	Count
40	149 x 165	+		90	55 x 18	+
41	150 x 37	+		91	59 x 116	+
42	152 x 42	+		92	6 x 131	+
43	153 x 108	+		93	61 x 162	+
44	153 x 20	?		94	63 x 151	not in box
45	154 x 144	-		95	7 x 46	-
46	156 x 166	+		96	8 x 61	not in box
47	16 x 4	+		97	Ler x Ler	+
48	16 x 66	+		98	Ler self	+
49	164 x 7	not in box		99	Cvi x Ler	+
50	166 x 25	+		100	Ler x Cvi	+